



Triana Science & Technology S.L.

Parque Tecnológico Ciencias de la Salud, Avenida Innovación, 1
18100 Armilla - Granada (Spain)
Tel: +34 958 750 583 - Fax: +34 958 750 584
Email: [triana@trianatech.com](mailto: triana@trianatech.com) WebSite: <http://www.trianatech.com>

CRYSTALLIZATION SCREEN FORMULATION

POLYETHYLEN GLYCOL 4000	Precipitant 1	[Precipitant 1]	Buffer	[Buffer]	pH
Box- 1	PEG 4000	30 %w/v	Sodium Acetate	100 mM	4.00
Box- 2	PEG 4000	30 %w/v	Sodium Acetate	100 mM	5.00
Box- 3	PEG 4000	30 %w/v	Tris-HCl	100 mM	6.00
Box- 4	PEG 4000	30 %w/v	Tris-HCl	100 mM	7.00
Box- 5	PEG 4000	30 %w/v	Tris-HCl	100 mM	8.00
Box- 6	PEG 4000	30 %w/v	Tris-HCl	100 mM	9.00